

# CD19-specific chimeric antigen receptor-modified T cells with safety switch produced under “Point-of-Care” using the *Sleeping Beauty* system for the very rapid manufacture and treatment of B-cell malignancies

Tim Chan<sup>1</sup>, John Gallagher<sup>1</sup>, Nai-Lin Cheng<sup>1</sup>, Fernando Carvajal-Borda<sup>1</sup>, Jacques Plummer<sup>1</sup>, Adeline Govekung<sup>1</sup>, John A. Barrett<sup>2</sup>, Pranay D. Khare<sup>2</sup>, Laurence J. N. Cooper<sup>2</sup>, Rutul R. Shah<sup>1</sup>

<sup>1</sup>Intrexon Corporation, Germantown MD; <sup>2</sup>ZIOPHARM Oncology, Boston, MA

## Abstract

Chimeric antigen receptor (CAR)-modified autologous T cells have shown remarkable clinical responses targeting B-cell malignancies. However, a majority of manufacturing approaches rely on expensive viral vectors and *ex vivo* activation of clinical-grade T cells for gene transfer with subsequent 2-3 week-long numeric expansion to achieve necessary dosing. We have engineered our clinically-validated *Sleeping Beauty* (SB) system for stable non-viral integration of CAR, cytokine signaling, and safety switch transgenes into T cells to eliminate need for viral vectors. We previously demonstrated that co-expression of membrane bound IL-15 (mIL15), Proc Natl Acad Sci U S A. 2016 Nov 29;113(48):E7788-E7797 on CD19-specific CAR<sup>+</sup> T cells using SB system allows for very rapid (<2 days) manufacturing with point-of-care (P-O-C) technology (Blood 2016 128:2807). This approach to production avoids the need for activation and propagation of CAR<sup>+</sup> T cells leading to augmented persistence and anti-tumor activity in an *in vivo* model of leukemia. We now update this technology by co-expressing a safety switch, HER11, on CD19-specific CAR<sup>+</sup> T cells. HER11 improves the safety profile by providing a mechanism for selective *in vivo* depletion of P-O-C-produced T cells through the administration of cetuximab, a clinically-available monoclonal antibody. We confirmed co-expression of CAR, mIL15, and HER11 in P-O-C-manufactured T cells and demonstrated that these cells exhibited a memory-like phenotype, potent antigen-specific cytotoxicity, and produced IFN $\gamma$ , TNF $\alpha$ , and IL-2 when co-cultured with CD19<sup>+</sup> tumor cells. Cetuximab specifically and efficiently killed mIL15<sup>hi</sup>HER11<sup>+</sup>CAR<sup>+</sup> T cells via antibody dependent cellular cytotoxicity and complement dependent cytotoxicity in a dose-dependent manner. Adoptive transfer of mIL15<sup>hi</sup>HER11<sup>+</sup>CAR<sup>+</sup> T cell generated in less than 2 days under P-O-C into immunocompromised (NSG) mice bearing established CD19<sup>+</sup> leukemia (NALM-6) resulted in significant reduction in tumor burden and improvement in overall survival, compared to control mice. Blood analysis demonstrated improved expansion and persistence of mIL15<sup>hi</sup>HER11<sup>+</sup>CAR<sup>+</sup> T cells compared to mIL15<sup>hi</sup>HER11<sup>-</sup>CAR<sup>+</sup> T cell control. In summary, we have generated T cells co-expressing mIL15 and the HER11 safety switch in less than two days using P-O-C that demonstrated superior efficacy and persistence compared to P-O-C T cells that expressed CAR without mIL15. These pre-clinical data demonstrate the effectiveness of very rapidly manufactured mIL15<sup>hi</sup>HER11<sup>+</sup>CAR<sup>+</sup> T cells in targeting CD19<sup>+</sup> tumors and provide a rationale for clinical evaluation.

## “Point-of-Care” (P-O-C) Manufacturing




Figure 1: Overview of the “Point-of-Care” (P-O-C) using the *Sleeping Beauty* (SB) system to very rapidly generate CAR<sup>+</sup> T cells. Very rapid manufacturing process removes need for *ex vivo* culture of T cells prior to infusion. Co-expression of mIL15 enhances *in vivo* persistence and anti-tumor response.

## Co-expression of CD19 CAR, mIL15 And Safety Switch In T Cells





Figure 2: Characterization of CAR+mIL15+HER11+ T cells. (A) Schematic of CD19-specific CAR and HER11 safety switch with mIL15 on cell surface of T cells genetically modified using DNA plasmids from the SB system. CD19-specific CAR is a chimeric fusion protein of an anti-human CD19 scFv and endodomains of chimeric CD28 and CD3 $\zeta$ . HER11 safety switch lacks signaling capability, but allows for conditional elimination of mIL15<sup>hi</sup> T cells using cetuximab via antibody dependent cellular cytotoxicity (ADCC). (B) Flow cytometric assessment of CAR, mIL15, and HER11 expression. For the purposes of *in vitro* assessments, T cells were repeatedly stimulated with feeder cells in presence of cytokines for 14 days to stably express these transgenes. Live cell population was gated based on FSC/SSC/CD45<sup>+</sup>/CD3<sup>+</sup>. (C) Confirmation of expression by Western blotting was confirmed by detection of chimeric CD3 $\zeta$  (left panel); HER11 was confirmed by co-immunoprecipitation followed by western blotting using anti-CD3 $\zeta$  antibodies (middle panel), and mIL15 expression was confirmed with anti-human IL-15 antibody (right panel).

## Implementing *Sleeping Beauty*-Modified Cells, a Non-viral Approach



## Antigen Specific Cytotoxicity And Cytokine Production By CD19 CAR+ T Cells




Figure 3: Specific cytotoxicity of CD19 CAR+mIL15+HER11+ T cells to CD19<sup>+</sup> target cell line. (A) Cytotoxicity measured by europium release assay using CD19<sup>+</sup> K562 (K562/CD19) and CD19<sup>-</sup> K562 (K562) cells at varying effector:target (E:T) ratios. (B) IFN $\gamma$ , TNF $\alpha$ , and IL-2 production following co-culture of CD19 CAR<sup>+</sup> T cells at different E:T ratios with the CD19<sup>+</sup> and CD19<sup>-</sup> target cells. Shown is the mean  $\pm$  SD from triplicate wells from a representative donor.

## Conditional Elimination of HER11 Expressing CD19 CAR+ T Cells




Figure 4: Ability of cetuximab to deplete CD19 CAR+mIL15+HER11+ T cells. (A) Schematic of ADCC. (B) Labeled CAR<sup>+</sup> T cells were treated with cetuximab or rituximab (control) for 24 hrs and co-cultured with allogeneic NK cells expressing Fc receptors. CAR<sup>+</sup> mIL15<sup>hi</sup> T cells remaining in culture post antibody treatment were determined by flow cytometry to calculate percent killing. Similar results were obtained using complement dependent cytotoxicity assays (data not shown). \*\*\* p<0.001; Student's T test.

## “Point-of-Care” CD19 CAR+ T Cells Eliminate Leukemia And Improve Survival In Mice




Figure 5: Anti-tumor activity of “Point-of-Care” CD19 CAR<sup>+</sup> T cells in CD19<sup>+</sup> B-cell (NALM-6) tumor model. (A) Study design. Immunocompromised (NSG) mice bearing established disseminated CD19<sup>+</sup> leukemias were treated on day 8 after tumor administration with saline, 10<sup>6</sup> mock transfected T cells, P-O-C CAR+mIL15<sup>hi</sup>HER11<sup>+</sup>, or P-O-C CAR+mIL15<sup>hi</sup>HER11<sup>+</sup> T cells. (B) Tumor burden serially measured by bioluminescence imaging (BLI) from the ventral view with the same 4 mice from respective groups. (C) Kaplan-Meier survival curves for overall survival. Arrows indicate the day genetically modified cells were administered. (D) Kaplan-Meier curves for overall survival. Arrows indicate the day genetically modified cells were administered. \*\*\* p<0.05; Log-Rank Test.

## P-O-C CD19 CAR+ T Cells Expand In Vivo And Persist At High Frequencies




Figure 6: *In vivo* expansion and persistence of “Point-of-Care” CD19 CAR<sup>+</sup> T cells in tumor bearing mice. (A) Representative flow cytometry data shown from whole blood evaluating expression of CAR and HER11 on day 22. Cell populations are gated based upon FSC/SSC/CD45<sup>+</sup>/CD3<sup>+</sup>. (B) Quantification of the CAR<sup>+</sup> T cells detected in whole blood. Shown is the mean  $\pm$  SEM, n = 4 to 11 mice/group.

## Summary

- T cells expressing CD19-specific CAR, mIL15 and safety switch were generated via very rapid manufacturing process (“Point-of-Care”) using the *Sleeping Beauty* system.
- “Point-of-Care” CAR<sup>+</sup> T cells reduces time to treatment by eliminating the need for *ex vivo* (i) activation and (ii) propagation of T cells.
- Safety switch provides a mechanism for conditional elimination post infusion.
- Co-expression of mIL15 enhances *in vivo* expansion and persistence of “Point-of-Care” CAR<sup>+</sup> T cells and improves anti-tumor response.
- “Point-of-Care” CD19 CAR<sup>+</sup> T cells eliminated CD19<sup>+</sup> tumor and improved survival in mice.
- These data support clinical evaluation of CAR<sup>+</sup> T cells co-expressing mIL15 and safety switch rapidly manufactured using the *Sleeping Beauty* system.

Conflict of Interest Disclosures: Chan, Intrexon Corporation, Employment; Gallagher, Intrexon Corporation, Employment; Cheng, Intrexon Corporation, Employment; Carvajal-Borda, Intrexon Corporation, Employment; Plummer, Intrexon Corporation, Employment; Govekung, Intrexon Corporation, Employment; Barrett, ZIOPHARM Oncology, Employment; Khare, ZIOPHARM Oncology, Employment; Cooper, Intrexon Corporation, Employment; Shah, Intrexon Corporation, Employment; Plummer, Intrexon Corporation, Employment; Khare, ZIOPHARM Oncology, Employment; Cooper, Intrexon Corporation, Employment; Shah, Intrexon Corporation, Employment; Plummer, Intrexon Corporation, Employment; Khare, ZIOPHARM Oncology, Employment; Cooper, Intrexon Corporation, Employment; Shah, Intrexon Corporation, Employment.